

## • Examination/Treatment Table 409

A Versatile Tool for Modern Nursing Environment



- **sturdy, functional construction**
- **versatile mattress base adjustments**
- **easy-to-clean surfaces**
- **wide selection of accessories**

The versatility of the Examination Table 409 by Merivaara makes it ideal for use during general examinations and medical and nursing procedures. The table features ergonomic, easy-to-use adjustments and a sturdy construction. The castors allow easy manouvering and lock securely in place. A wide selection of accessories add to the functions of the table during nursing procedures.

# 409 Examination/Treatment Table

## Technical Specifications

- Frame of epoxy-powder-coated steel
- Mattress base 2-sectional, upholstery of artificial leather, easily cleanable
- Sleepless back section adjustment with gas spring and height adjustment hydraulically (590-260 mm) or height and back section adjustments electrically
- Height adjustment pedals on either side of the table
- Castors Ø 125 mm with brakes, locks directionally
- Push bars on both ends

## Dimensions

- Max. length 2030 mm
- Max. width 705 mm

## Accessories and Options

- paper roll holder 732
- paper cutting mechanism 734
- wire basket 583
- trendelenburg and antitrendelenburg adjustments 977
- castors with central braking system 6321



Examination table 409 with trendelenburg adjustment.



Examination table 409 with available accessories: paper roll holder and cutting mechanism, wire basket, central braking system and trendelenburg adjustments.